



## CONTACT US

[www.vathin.com](http://www.vathin.com)  
[marketing@vathin.com](mailto:marketing@vathin.com)

## Single-use Flexible Ureteroscope

Urology solution



# Single-use Flexible Ureteroscope



Single-use Flexible Ureteroscope system is specially designed to meet the evolving needs of urology specialists, providing a comprehensive all-in-one-solution for healthcare professionals and offering a range of advantages that make it the ideal choice for clinical application. With its unique features, the system ensures smoother drainage and clearer picture quality, enabling urologists to perform procedures with utmost accuracy.

## Main Features >>




### Suction

Smoother drainage, clearer picture  
Easier to operate with rotatable suction tube.

### Bending lock

More convenient operation, more accurate positioning



### 90° Horizontal Rotation


Optimize the view without grip adjustment

### Riveted Snake bone


Providing flexibility & durability for better performance  
Bending capability  
Ureteroscope: 0~285° Cystoscope: 0~210°



## Specifications




Model	Working length (mm)	Working channel ID (mm/Fr)	Outer diameter of main hose OD (mm/Fr)	Bending section	Field of view	Orientation	Bending lock	Suction function	Color
US-S170	700	1.2/3.6Fr	2.8/8.4Fr	U285°/D285°	110°±10%	Standard	Y	Y	Grey
US-S170R		1.2/3.6Fr	2.8/8.4Fr			Standard			
US-S180		1.4/4.2Fr	3.1/9.3Fr			Reverse			
US-E170		1.2/3.6Fr	2.8/8.4Fr			Standard	N	N	Grey
US-E170R		1.2/3.6Fr	2.8/8.4Fr			Standard			
US-E180		1.4/4.2Fr	3.1/9.3Fr			Reverse			



### Reverse deflection:

Moving the lever down bends the tip up; moving it up bends the tip down.



### Standard deflection:

Moving the lever down bends the tip down; moving it up bends the tip up.